


APPENDIX/INDEX

Selecting and Designing a Filtration System	120
Graphs: Flow Rate vs Pressure	121
Chemical Compatibility	122
Membrane Filters	122
Disposable Syringe Filter Units	124
Capsules	126
Cartridges	128
Cartridges and O-rings	130
Sterile Membrane Cross Reference Conversion Guide	131
Glass, Quartz, and Paper Cross Reference Conversion Guide	131
Index by Product Name	132
Index by Catalog Number	135



Selecting and Designing a Filtration System

In designing or selecting a system for microporous filtration, it is worthwhile to spend some time defining a few specific parameters of the filtration problem:

- **What is to be filtered, liquid or gas?** For filtering of gases or sterile venting of vessels or gas lines, use a hydrophobic membrane such as PTFE. Chemical resistance tables can be used to help select a membrane material suitable to the liquid being filtered. Also check the compatibility of the filter holder and any O-rings with the filtrate.
- **What size are the smallest particles to be removed?** In general, a pore size just smaller than the smallest size particle to be removed is chosen. For example, 0.20 μm pore size membrane is generally used for sterilizing culture medium, but a 0.10 μm membrane is needed to ensure Mycoplasma removal.
- **At what pressures will the system be operated?** Check the specifications or assume a differential pressure of 2-5 psi ($\sim 0.3 \text{ kg/cm}^2$) for a single holder. Be sure to account for the effects of pressure drop in a multi-holder system.
- **How will the filter holder be connected to the system?** Standard Luer, threaded, sanitary, and hose connectors are available depending on the model. Check the specifications for the product you choose.
- **What is the operating temperature?** Check the specifications for your filtration medium (membrane or other filter) and holder.
- **What size filter and holder do I need?** Estimate what volume you will be filtering and the amount of time you plan to allow for the filtration. By dividing the volume by the time in which you expect to complete this filtration you will get a desired minimum flow rate. Use the graphs at right to determine the flow rate per unit area for the differential pressure and membrane in your system (assume a pressure differential of 5 psi ($\sim 0.3 \text{ kg/cm}^2$) if it is not known). Divide the desired flow rate by the flow rate per unit area to get a minimum filtration area in cm^2 . Multiply this area by a factor of ~ 5 to allow for clogging.

ADJUSTMENTS

For gaseous filtrations that take place outside of standard temperature and pressure (20°C , one atmosphere), correct flow rates using the following formula:

$$F = F_o \frac{293}{273+T} \frac{P + \Delta P/2}{1.03 + \Delta P/2}$$

F = Corrected flow rate
 F_o = Flow rate from chart
 T = Temperature in $^\circ\text{C}$
 P = Exit pressure in kg/cm^2
 ΔP = Pressure drop through the system in kg/cm^2


Liquid viscosity can have a significant effect on flow rate. Use this formula to correct for this effect:

$$\text{Flow rate} = \frac{A \cdot P}{V}$$


A = Effective filtration area
 P = Differential pressure
 V = Viscosity

GRAPHS – FLOW RATE VS. DIFFERENTIAL PRESSURE


MIXED ESTER OF CELLULOSE (AIR)




MIXED ESTER OF CELLULOSE (WATER)



CELLULOSE ACETATE (AIR)



CELLULOSE ACETATE (WATER)



Chemical Compatibility

This guide should be used as a reference to select the proper membrane, O-ring and filter holder with regard to chemical compatibility.

Recommendations are based on a 24-hour static exposure to the test fluid at room temperature.

- O – Recommended.** No change in either water flow rate or bubble point will be observed.
- * – Limited resistance.** Additional in-house testing is advised as swelling, discoloration or other minor changes may occur.
- x – Not recommended.** Significant changes in water flow rate and/or bubble point can be expected.
- – Data not available.**

CHEMICAL COMPATIBILITY – MEMBRANE FILTERS

Chemical		Polymer/Product											
		MCE		CA	NYLON	PCTE		PTFE	Supported PTFE	Hydrophilic PTFE	CMF Coated Cellulose acetate	37 mm monitor	37 mm monitor set
		White	Black			White	Black						
ACID	Glacial Acetic acid	X	X	X	X	O	X	O	O	O	X	X	X
	10% Acetic acid	O	O	O	O	O	O	O	O	O	O	*	*
	12 kmol/m ³ Hydrochloric acid (37%, 12N)	X	X	X	X	O	X	O	O	O	X	X	X
	6 kmol/m ³ Hydrochloric acid (19%, 6N)	O	X	X	X	O	O	O	O	O	X	*	*
	12 kmol/m ³ Nitric acid (53%, 12N)	X	X	X	X	O	X	O	O	O	X	X	X
	6 kmol/m ³ Nitric acid (26%, 6N)	O	X	X	X	O	X	O	O	O	X	X	X
	18 kmol/m ³ Sulfuric acid (96%, 36N)	X	X	X	X	X	X	O	O	X	X	X	X
	3 kmol/m ³ Sulfuric acid (16%, 6N)	O	O	X	X	O	O	O	O	O	X	*	*
	85% Phosphoric acid	O	O	X	-	O	X	O	O	O	X	*	*
	5% Boric acid	O	O	O	*	O	O	O	O	O	O	O	O
	50% Formic acid	*	X	X	-	O	X	O	O	O	*	-	-
35% Hydrofluoric acid	O	X	X	X	O	X	O	O	O	X	X	X	
60% Perchloric acid	X	X	X	-	O	X	O	O	O	X	-	-	
ALKALI	6 kmol/m ³ Sodium hydroxide (26%, 6N)	X	X	X	X	X	X	O	O	O	X	X	X
	6 kmol/m ³ Potassium hydroxide (20%, 6N)	X	X	X	O	X	X	O	O	O	X	X	X
	6 kmol/m ³ Aqueous ammonia (11%, 6N)	X	X	X	X	X	X	O	O	O	X	X	X
ALCOHOL	Methyl alcohol	X	X	*	O	O	X	O	O	O	O	X	X
	Ethyl alcohol	X	X	*	O	O	X	O	O	O	O	X	X
	Isopropyl alcohol	*	X	*	O	O	X	O	O	O	O	X	X
	Isobutyl alcohol	O	X	O	O	O	X	O	O	O	O	X	X
	Butyl alcohol	*	X	O	O	O	X	O	O	O	O	X	X
	Glycerol	O	O	O	O	O	O	O	O	O	O	O	O
	Amyl alcohol	*	X	O	O	O	X	O	O	O	O	*	*
	Benzyl alcohol	X	X	O	O	O	X	O	O	O	X	*	*
Ethylene glycol	*	X	O	O	O	O	O	O	O	O	*	X	

Continued on next page

CHEMICAL COMPATIBILITY – MEMBRANE FILTERS (CONTINUED)

Chemical		Polymer/Product											
		MCE		CA	NYLON	PCTE		PTFE	Supported PTFE	Hydrophilic PTFE	CMF Coated Cellulose acetate	37 mm monitor	37 mm monitor set
		White	Black			White	Black						
ETHERS	Ethyl ether	○	○	○	○	*	x	○	○	○	○	*	x
	Isopropyl ether	○	○	○	-	○	x	○	*	○	○	x	x
	Tetrahydrofuran	x	x	x	*	x	x	○	*	○	x	x	x
	Dioxane	x	x	x	○	x	x	○	*	○	x	x	x
	Petroleum ether	○	○	○	-	○	○	○	○	○	○	-	-
ESTERS	Methyl acetate	x	x	x	○	x	x	○	*	○	x	x	x
	Butyl acetate	x	x	*	-	*	x	○	*	○	○	x	x
	Amyl acetate	x	x	*	○	○	*	○	*	○	○	x	x
KETONES	Acetone	x	x	x	○	x	x	○	○	○	x	x	x
	Methylethyl ketone	x	x	x	○	x	x	○	*	○	x	x	x
	Methyl isobutyl ketone	x	x	x	-	x	x	○	*	○	x	x	x
	Cyclohexanone	x	x	x	○	x	x	○	○	○	x	x	x
HYDROCARBONS	Benzene	○	x	○	*	x	x	○	○	○	○	x	x
	Toluene	○	○	○	*	x	*	○	*	○	○	x	x
	Xylene	○	x	○	○	○	x	○	*	○	x	x	x
	Hexane	○	○	○	○	○	○	○	○	○	○	x	x
	Gasoline	○	○	○	○	○	○	○	*	○	○	*	x
	Kerosene	○	○	○	○	○	○	○	*	○	x	○	*
HALOGENATED HYDROCARBONS	Chloroform	○	○	x	○	x	x	○	*	○	x	x	x
	Methylene chloride	*	x	x	*	x	x	○	*	○	x	x	x
	Trichloroethylene	○	○	○	○	x	x	○	*	○	○	x	x
	Tetrachloroethylene	○	x	○	-	x	x	○	*	○	○	x	x
	Carbon tetrachloride	○	○	○	○	○	○	○	○	○	○	x	x
AMINES	Aniline	x	x	x	x	x	x	○	○	○	x	x	x
	Dimethyl formamide	x	x	x	○	x	x	○	○	○	x	x	x
	Diethyl acetamide	x	x	x	○	x	x	○	○	○	x	x	x
	Triethanolamine	x	x	○	○	○	x	○	○	○	x	x	x
CELLOSOLVES	Methyl cellosolve	x	x	x	○	○	x	○	○	○	x	x	x
	Butyl cellosolve	x	x	○	○	○	x	○	○	○	○	x	x
MISCELLANEOUS	Nitrogen	○	○	○	○	○	○	○	○	○	○	○	○
	Hydrogen	○	○	○	○	○	○	○	○	○	○	○	○
	30% Hydrogen peroxide	○	○	○	○	○	○	○	○	○	○	x	x
	Saline solution	○	○	○	○	○	○	○	○	○	○	○	○
	Dimethylsulfoxide	x	x	*	*	x	x	○	○	○	x	x	x
	Nitrobenzene	x	x	x	○	x	x	○	○	○	x	x	x
	Methanol (1): Chloroform (1)	x	x	x	-	x	x	○	*	○	x	x	x
	Pyridine	x	x	x	○	x	x	○	*	○	x	x	x
	Acetonitrile	x	x	x	○	x	x	○	○	○	x	x	x
	Phenol	x	x	x	○	x	x	○	○	○	x	x	x
	Freon	○	○	○	○	○	○	○	*	○	○	○	○
	37% Formaldehyde	○	x	*	○	○	x	○	○	○	○	x	x
	Silicone oil	○	○	○	○	○	x	○	*	○	x	-	-
	n-Hexane (95): Ethyl acetate (5)	*	○	x	-	○	○	○	*	○	○	x	x
	Nitric acid (70): Distilled water (30)	x	x	x	x	x	x	○	○	○	x	x	x
Petroleum oil	○	○	○	○	○	x	○	○	○	○	*	x	

○ – Recommended * – Limited resistance x – Not recommended - – Data not available

CHEMICAL COMPATIBILITY: DISPOSABLE SYRINGE FILTER UNITS

Chemical		AS Mixed cellulose esters with acrylic housing	CS Cellulose acetate with acrylic housing	CP Cellulose acetate with PP housing	HP Hydrophilic PTFE with PP housing	JP Hydrophobic PTFE with PP housing	NS Nylon with acrylic housing	NP Nylon with PP housing
ACIDS	3 kmol/m ³ Hydrochloric acid (10%, 3N)	O	*	*	O	O	*	*
	9 kmol/m ³ Hydrochloric acid (30%, 9N)	X	X	X	O	O	X	X
	1 kmol/m ³ Sulfuric acid (5%, 2N)	O	O	O	O	O	*	*
	4 kmol/m ³ Sulfuric acid (20%, 8N)	X	X	X	O	O	X	X
	1 kmol/m ³ Nitric acid (5%, 1N)	O	O	*	O	O	*	*
	5 kmol/m ³ Nitric acid (20%, 5N)	X	X	X	*	*	X	X
	20% Acetic acid	O	O	O	O	O	O	O
	Glacial acetic acid	X	X	X	O	O	X	X
	10% Hydrofluoric acid	X	X	X	O	O	X	X
	35% Hydrofluoric acid	X	X	X	O	O	X	X
	10% Chromic acid	*	*	*	*	*	*	*
10% Phosphoric acid	O	O	O	O	O	*	*	
ALKALIS	2.5 kmol/m ³ Sodium hydroxide (10%, 2.5N)	X	X	X	O	O	X	X
	2 kmol/m ³ Potassium hydroxide (10%, 2N)	X	X	X	O	O	X	O
	8 kmol/m ³ Aqueous ammonia (28%, 8N)	X	O	*	O	O	O	O
ALCOHOLS	Methyl alcohol	X	X	O	O	O	O	O
	Ethyl alcohol	X	*	O	O	O	O	O
	n-Propyl alcohol	*	O	O	O	O	O	O
	Isopropyl alcohol	*	O	O	O	O	O	O
	n-Butyl alcohol	*	O	O	O	O	O	O
	Amyl alcohol	*	O	O	O	O	O	O
	Benzyl alcohol	*	X	X	O	O	*	O
	Ethylene glycol	X	X	O	O	O	*	O
Glycerol	O	O	O	O	O	O	O	
ETHERS	Ethyl ether	X	X	*	*	*	X	O
	Isopropyl ether	X	X	O	O	O	-	-
	Tetrahydrofuran (THF)	X	X	X	*	*	X	*
	Dioxane	X	X	X	*	*	X	O
ESTERS	Methyl acetate	X	X	X	*	*	X	O
	Ethyl acetate	X	X	X	*	*	X	-
	Butyl acetate	X	X	X	*	*	X	X
	Amyl acetate	X	X	*	*	*	X	-
KETONES	Acetone	X	X	X	O	O	X	O
	Methyl ethyl ketone (MEK)	X	X	X	*	*	X	O
	Methyl isobutyl ketone (MIBK)	X	X	X	*	*	X	O
	Cyclohexanone	X	X	X	*	*	X	O
HYDROCARBONS	Benzene	X	X	*	*	*	X	*
	Toluene	X	X	*	*	*	X	*
	Xylene	X	X	*	*	*	X	*
	n-Hexane	X	X	*	*	*	X	O
	Gasoline	*	*	O	O	O	*	O
	Kerosene	O	O	*	*	*	X	O

Continued on next page

CHEMICAL COMPATIBILITY: DISPOSABLE SYRINGE FILTER UNITS (CONTINUED)

Chemical		AS Mixed cellulose esters with acrylic housing	CS Cellulose acetate with acrylic housing	CP Cellulose acetate with PP housing	HP Hydrophilic PTFE with PP housing	JP Hydrophobic PTFE with PP housing	NS Nylon with acrylic housing	NP Nylon with PP housing
HALOGENATED HYDROCARBONS	Chloroform	X	X	X	*	*	X	*
	Methylene chloride	X	X	X	*	*	X	*
	Trichloroethylene	X	X	*	*	*	X	O
	Carbon tetrachloride	X	X	*	O	O	X	O
	Trichloroethane	X	X	X	*	*	X	O
	Perchloroethylene	X	X	X	*	*	X	O
	Freon (TMC)	X	X	*	*	*	X	O
AMINES	Aniline	X	X	X	*	*	X	*
	Dimethyl formamide	X	X	X	*	*	X	O
	Diethylacetamide	X	X	X	*	*	X	-
	Triethanolamine	X	O	O	O	O	X	-
CELLOSOLVES	Ethyl acetate cellosolve	X	X	X	*	*	X	*
MISCELLANEOUS	Acetonitrile	X	X	X	*	*	X	*
	Pyridine	X	X	X	*	*	X	O
	6% Sodium Hypochloride	X	X	X	O	O	X	-
	35% Formaldehyde	X	*	*	O	O	X	-
	Iron (II) chloride	O	O	O	O	O	O	O
	Copper sulfate	O	O	O	O	O	O	O
	Mineral oil	*	O	*	*	*	*	*
	Salt water	O	O	O	O	O	O	O
	10% Hydrogen peroxide	X	*	O	O	O	X	O
	Nitrobenzene	X	X	X	*	*	X	-
	Phenol	X	X	X	O	O	X	-
	Silicone oil	X	X	X	O	O	X	-
	Petroleum oil	*	*	O	O	O	X	-
	Acetonitrile (70): water (30)	X	X	X	O	O	X	O

O – Recommended

* – Limited resistance

X – Not recommended

- – Data not available

CHEMICAL COMPATIBILITY – CAPSULES

Chemical	CCS	CCF	CCFH	CCP	CCG	
ACIDS	5% Acetic acid	O	O	O	O	O
	20% Acetic acid	O	O	*	O	*
	10% Chromic acid	-	*	*	*	X
	3 kmol/m ³ Hydrochloric acid (10%, 3N)	O	O	O	O	O
	11 kmol/m ³ Hydrochloric acid (35%, 11N)	O	O	X	O	*
	10% Hydrofluoric acid	-	O	O	O	X
	1 kmol/m ³ Nitric acid (5%, 1N)	-	O	O	O	*
	4 kmol/m ³ Nitric acid (20%, 4N)	-	*	*	*	X
	10% Phosphoric acid	-	O	O	O	O
	1 kmol/m ³ Sulfuric acid (5%, 2N)	O	O	O	O	O
4 kmol/m ³ Sulfuric acid (20%, 8N)	O	O	O	O	*	
ALKALIS	2.5 kmol/m ³ Sodium hydroxide (10%, 3N)	O	O	X	O	*
	2 kmol/m ³ Potassium hydroxide (10%, 2N)	O	O	X	O	*
	5 kmol/m ³ Aqueous ammonia (10%, 5N)	O	O	X	O	*
	15 kmol/m ³ Aqueous ammonia (28%, 15N)	O	O	X	O	*
ALCOHOLS	Methyl alcohol	O	O	O	O	*
	Ethyl alcohol	O	O	O	O	*
	n-propyl alcohol	O	O	O	O	*
	n-butyl alcohol	O	O	*	O	*
	Ethylene glycol	O	O	O	O	O
ETHERS	Ethyl ether	*	*	*	*	X
	Dioxane	*	*	*	*	X
	Tetrahydrofuran (THF)	X	*	X	*	X
ESTERS	Amyl acetate	-	*	*	*	-
	Methyl acetate	-	*	*	*	*
	Ethyl acetate	X	*	*	*	*
	Butyl acetate	-	*	*	*	*
KETONES	Acetone	X	*	*	*	X
	Methyl ethyl ketone (MEK)	X	*	*	*	X
	Methyl isobutyl ketone (MIBK)	X	*	-	*	X
HYDROCARBONS	n-hexane	*	*	-	*	*
	Cyclohexane	*	*	*	*	*
	Benzene	*	*	*	*	*
	Toluene	X	*	*	*	X
	Xylene	-	*	*	*	*
HALOGENATED HYDROCARBONS	Chloroform	X	*	X	*	X
	Carbon tetrachloride	-	*	*	*	*
	Freon (TMC)	X	*	-	*	X
	Methylene chloride	X	*	X	*	X
	Trichloroethylene	-	*	*	*	X
	Trichloroethane	-	*	-	*	*
AMINES	Dimethyl formamide	X	*	*	*	X
ALDEHYDES	Acetaldehyde	-	-	-	-	-
	35% Formaldehyde	O	O	*	O	O

Continued on next page

CHEMICAL COMPATIBILITY – CAPSULES (CONTINUED)

Chemical		CCS	CCF	CCFH	CCP	CCG
CELLOSOLVES	Acetic cellosolve	-	*	-	*	*
MISCELLANEOUS	Acetonitrile	-	*	-	*	-
	Pyridine	-	*	-	*	X
	Nitrobenzene	-	-	X	-	-
	6% Sodium hypochlorite	-	O	O	O	O
	Ferrous chloride	O	O	O	O	O
	Copper sulfate	O	O	O	O	O
	Mineral oil	*	*	*	*	*
	Salt water	O	O	O	O	O
	10% Hydrogen peroxide	O	O	O	O	*

O – Recommended * – Limited resistance X – Not recommended - – Data not available

CHEMICAL COMPATIBILITY – CARTRIDGES

Chemicals		TCR	TCS TCS-G TCS-E	TCF TCFH	TCY TCYE	TCP TCPE TCPD	TC	TC (SUS type with NBR Gasket)	TCG	TCG-R (NBR Gasket)
ACIDS	20% Acetic acid	O	O	O	O	O	X	*	*	*
	10% Chromic acid	X	-	*	X	*	X	X	X	X
	3 kmol/m ³ Hydrochloric acid (10%, 3N)	X	O	O	X	O	X	X	O	X
	9 kmol/m ³ Hydrochloric acid (30%, 9N)	X	O	O	X	O	X	X	*	X
	10% Hydrofluoric acid	X	-	O	X	O	*	X	X	X
	1 kmol/m ³ Nitric acid (5%, 1N)	X	-	O	*	O	*	X	*	X
	4 kmol/m ³ Nitric acid (20%, 4N)	X	-	*	X	*	X	X	X	X
	10% Phosphoric acid	O	-	O	O	O	*	X	O	O
	1 kmol/m ³ Sulfuric acid (5%, 2N)	X	O	O	O	O	*	*	O	*
4 kmol/m ³ Sulfuric acid (20%, 8N)	X	O	O	X	O	X	X	*	X	
ALKALIS	2.5 kmol/m ³ Sodium hydroxide (10%, 2.5N)	X	O	O	X	*	X	X	X	X
	2 kmol/m ³ Potassium hydroxide (20%, 2N)	X	O	O	X	*	X	X	X	X
	15 kmol/m ³ Aqueous ammonia (28%, 15N)	*	O	O	X	*	X	X	*	*
ALCOHOLS	Methyl alcohol	*	O	O	O	O	*	*	*	*
	Ethyl alcohol	O	O	O	O	O	*	*	*	*
	n-propyl alcohol	O	O	O	O	O	O	*	*	*
	n-butyl alcohol	O	O	O	O	O	O	*	*	*
ETHERS	Dioxane	X	X	*	X	*	X	X	X	X
	Tetrahydrofuran (THF)	X	X	*	X	*	X	X	X	X
	Ethyl ether	X	X	X	X	X	X	X	X	X
ESTERS	Amyl acetate	-	-	X	X	X	-	-	-	-
	Methyl acetate	X	-	*	X	*	*	X	*	X
	Ethyl acetate	X	X	*	X	*	*	X	*	X
	Butyl acetate	X	-	*	X	*	*	X	*	X
KETONES	Acetone	X	X	*	X	*	-	X	X	X
	Methyl ethyl ketone (MEK)	X	X	*	X	*	X	X	X	X
	Methyl isobutyl ketone (MIBK)	X	X	*	X	*	X	X	X	X
HYDROCARBONS	n-hexane	X	X	X	X	X	X	X	X	X
	Benzene	X	X	*	*	*	X	X	*	X
	Toluene	X	X	*	*	*	*	X	X	X
	Xylene	X	-	*	*	*	*	X	*	X
HALOGENATED HYDROCARBONS	Chloroform	X	X	X	X	X	X	X	X	X
	Carbon tetrachloride	*	-	*	*	*	*	X	*	X
	Freon (TMC)	X	X	*	X	*	-	-	X	X
	Methylene chloride	X	X	X	X	X	X	X	X	X
	Trichloroethylene	X	-	X	X	X	X	X	X	X
	Trichloroethane	X	-	*	*	*	*	X	*	-
AMINES	Dimethyl formamide	X	X	*	X	*	X	X	X	X

Continued on next page

CHEMICAL COMPATIBILITY – CARTRIDGES (CONTINUED)

Chemicals		TCR	TCS TCS-G TCS-E	TCF TCFH	TCY TCYE	TCP TCPE TCPD	TC	TC (SUS type with NBR Gasket)	TCG	TCG-R (NBR Gasket)
ALDEHYDES	35% Formaldehyde	*	O	O	*	O	*	*	O	*
CELLOSOLVES	Acetic cellosolve	X	-	*	X	*	-	-	*	-
MISCELLANEOUS	Acetonitrile	X	-	*	X	*	-	-	-	-
	Pyridine	X	-	*	X	*	X	X	X	X
	6% Sodium hypochlorite	X	-	*	X	*	X	X	*	*
	Ferrous chloride	O	O	O	O	O	X	X	O	X
	Copper sulfate	O	O	O	O	O	*	O	O	O
	Mineral oil	*	*	*	*	X	X	X	*	*
	Salt water	O	O	O	O	O	O	*	O	*
	10% Hydrogen peroxide	*	O	O	*	O	*	*	*	*

O – Recommended * – Limited resistance X – Not recommended - – Data not available

CHEMICAL COMPATIBILITY – CARTRIDGES AND O-RINGS

Chemicals		Gasket/O-ring										
		WPP	WPS	WCP	WSP	EPR EPDM	SILICONE	NBR	CHLOROPRENE	FPM	PTFE	FEP Encap-sulated
ACIDS	20% Acetic acid	○	*	*	*	○	○	*	*	○	○	○
	10% Chromic acid	*	*	X	○	*	*	X	X	*	○	○
	3 kmol/m ³ Hydrochloric acid (10%, 3N)	○	X	X	X	○	○	○	○	○	○	○
	9 kmol/m ³ Hydrochloric acid (30%, 9N)	○	X	X	X	○	X	*	X	○	○	○
	10% Hydrofluoric acid	○	X	X	X	*	X	X	*	*	○	○
	1 kmol/m ³ Nitric acid (5%, 1N)	○	*	*	*	○	○	X	X	○	○	○
	4 kmol/m ³ Nitric acid (20%, 4N)	*	*	X	*	X	X	X	X	○	○	○
	10% Phosphoric acid	○	○	*	○	○	○	○	○	○	○	○
	1 kmol/m ³ Sulfuric acid (5%, 2N)	○	*	*	*	○	○	○	○	○	○	○
4 kmol/m ³ Sulfuric acid (20%, 8N)	○	X	X	X	○	○	○	○	○	○	○	
ALKALIS	2.5 kmol/m ³ Sodium hydroxide (10%, 2.5N)	○	○	*	○	○	X	○	○	○	○	○
	2 kmol/m ³ Potassium hydroxide (20%, 2N)	○	○	*	-	○	*	○	○	○	○	○
	15 kmol/m ³ Aqueous ammonia (28%, 15N)	○	○	*	*	○	○	X	○	○	○	○
ALCOHOLS	Methyl alcohol	○	○	○	○	○	○	X	○	*	○	○
	Ethyl alcohol	○	○	○	○	○	○	*	○	○	○	○
	n-propyl alcohol	○	○	○	-	○	○	*	○	○	○	○
	n-butyl alcohol	○	○	○	*	○	○	X	○	○	○	○
ETHERS	Dioxane	*	*	-	-	*	*	X	X	X	○	○
	Tetrahydrofuran (THF)	*	*	X	-	X	X	X	X	X	○	○
	Ethyl ether	*	*	○	-	X	X	X	*	X	○	○
ESTERS	Amyl acetate	*	*	-	-	*	X	X	X	X	○	○
	Methyl acetate	*	*	○	-	○	*	X	*	X	○	○
	Ethyl acetate	*	*	○	*	*	*	X	X	X	○	○
	Butyl acetate	*	*	○	-	*	X	X	X	X	○	○
KETONES	Acetone	*	*	*	○	○	*	X	*	X	○	○
	Methyl ethyl ketone (MEK)	*	*	○	*	*	*	X	X	X	○	○
	Methyl isobutyl ketone (MIBK)	*	*	*	-	X	○	X	X	X	○	○
HYDROCARBONS	n-hexane	*	*	○	-	X	X	○	*	○	○	○
	Benzene	*	*	*	*	X	*	X	X	○	○	○
	Toluene	*	*	○	*	X	X	X	X	○	*	*
	Xylene	*	*	○	○	X	X	X	X	○	○	○
HALOGENATED HYDROCARBONS	Chloroform	*	*	○	*	X	X	X	X	○	*	*
	Carbon tetrachloride	*	*	○	-	X	X	X	X	○	○	○
	Freon (TMC)	*	*	○	-	-	-	X	-	-	○	○
	Methylene chloride	*	*	○	-	X	X	X	X	○	○	○
	Trichloroethylene	*	*	○	*	X	X	X	X	○	○	○
Trichloroethane	*	*	○	-	X	X	X	-	○	○	○	
AMINES	Dimethyl formamide	X	*	*	*	○	○	X	X	X	○	○
ALDEHYDES	35% Formaldehyde	○	*	*	*	○	-	○	○	○	○	○
CELLOSOLVES	Acetic cellosolve	*	*	-	-	*	○	X	*	X	○	○
MISCELLANEOUS	Acetonitrile	*	*	○	-	○	X	X	○	X	○	○
	Pyridine	*	-	○	-	X	X	X	X	X	○	○
	6% Sodium hypochlorite	*	X	X	X	○	X	*	○	○	○	○
	Ferrous chloride	○	X	X	X	○	○	○	○	*	○	○
	Copper sulfate	○	○	○	○	○	○	○	○	○	○	○
	Mineral oil	*	*	○	○	X	*	○	*	○	○	○
	Salt water	○	*	*	*	○	○	○	○	○	○	○
10% Hydrogen peroxide	*	○	*	○	*	○	X	○	○	○	○	

Sterile Membrane Cross Reference Conversion Guide

Advantec MFS	Description	Millipore	Whatman	Pall	Catalog Page
A020H047A	White grid	GSWG 047 S1	7187 114	-	19
A045C047A	White plain	HAWP 047 S0	-	-	19
A045D047A	White grid	HAWG 047 S0	-	63077	19
A045F047A	White grid	HAWG 047 S2, SJ	7141 104	-	19
A045F047W	White grid	HAWG 047 S2	-	-	19
A045H047A	White grid	HAWG 047 S1, S	7141 114	66068, 66278	19
A045H047W	White grid	HAWG 047 S1, S	7141 114	66068, 66278	19
A045H047Y	White grid	HAWG 047 S3	7141 124	66278	19
A045T047A	White grid	HAWG 047 A0	7141 204	64194	19
A065F047A	White grid	HCWG 047 S4	-	-	19
A065H047A	White grid	HCWG 047 S1	-	66426	19
A045R047A	Black grid	HABG 047 S1, S	7153 104	66378	20
A080R047A	Black grid	AABG 047 SO	-	-	20
A045W047A	Green grid	-	7155 104	66379	20

See individual listings for packaging configurations.

Glass, Quartz, and Paper Cross Reference Conversion Guide

Advantec MFS	Millipore	Whatman	Pall	Fisher	Schleicher & Schuell	Ahlstrom	Catalog page
GA-55	APFA	GF/A	-	09-804-55A	31	111	28
GF-75	APFF	GF/F	-	-	20	151	28
GB-100R	-	EPM2000	(Air)	-	1HV	-	28
GB-140	APFB	GF/B	-	09-804-55B	32	121	28
GC-50	AP-40/APFC	GF/C 934 AH	A/E (Water)	09-804-55C	30/25	131	28
GC-90	AP15	-	-	-	-	-	28
GD-120	APFD	GF/D	-	09-804-55D	40	141	28
GS-25	AP20	-	-	-	-	164	28
QR-100	-	QM-A	Micro Quartz	-	-	-	28
No. 1	FP105	1	-	09-795F	410 or 1450cv	631	32
No. 101	FH560-260	4	-	09-790-12F	-	-	32
No. 131	FP102	3 or 6	-	09-803-6G	597	238	32
No. 2	FP104	2	-	09-801-E	604	601	32
No. 231	-	1	-	-	-	-	32
No. 232	-	-	-	-	-	642	32
No. 235	-	5	-	09-803-5F	-	610	32

Comparisons of glass fiber filters are for those with closest similarity. Exact counterparts do not exist due to the nature of the depth type filter matrix.

A

Absorbent pads	.21
Accessories, glassware	.86
Acetate membrane	.7, 20
Acid/alkali test papers	.44
Aerosol filter holder	.40, 102
Air filters	.39-41
Air Monitoring Cassettes	.40
All glass filtration apparatus	.87
Ampouled media	.22
Asbestos monitoring	.41
Ashless filter paper	.33
Autoclave pack filters	.19

B

Bacteriological filters	.18-20
Black filters	.4, 6, 14, 20
Blaine test papers	.44
Blood filter, Nobuto	.43
Blotting membrane	.35
Borosilicate glass filter holders	.79-88
Borosilicate glass filters	.26, 28-29

C

Capsule filters	.48-52
Glass	.52
PES	.49
Polypropylene	.51
PTFE	.50
Cartridge filters	.53-70
Cartridge housings	.71-76
Cellulose acetate	.55
Coated cellulose acetate	.61-62
Epoxy cellulose	.66
Glass fiber	.68
Melt blown polypropylene	.69
Pleated	.55-64, 66-68
Polyethersulfone (PES)	.56-58
Polypropylene	.63-65, 69
PTFE	.59-60
Spiral wound	.65
String wound	.70
Cassettes	.40-41
Cellulosic materials	
Cellulose acetate	.7, 15, 20, 55
Mixed cellulose esters (MCE)	.4-6, 15, 19-20, 40-41
Nitrocellulose, pure	.35
Prefilters	.25
Silicone treated	.38

Support pads	.41
Thimbles	.42
Characteristics of filters	
Depth	.24
Membrane	.2-3
Chemical compatibility guides	
37mm monitors	.122-123
Capsule filters	.126-127
Cartridge filters	.128-130
Disposable filter units	.124-125
Membrane filters	.122-123
O-rings	.130
Chromatography grade papers	.36
CMF	.9, 61-62
Coated membrane	.9, 61-62
Coliform analysis	.18
Composite filter	.39
Conversion guides	.131
Creped paper	.37
Culture media	.22

D

Depth filters	.24-30
Cellulose	.25
Glass fiber	.26, 28-29
PTFE	.30
Quartz	.27-29
Disposable filter units	.15-16, 45

E

Edge hydrophobic filters	.6, 40
Electrophoresis	.35-36
Epifluorescence microscopy	.14
Escherichia coli	.18
Extraction manifold	.93
Extraction thimbles	.42

F

Fecal coliform analysis	.18
Fibrous filters	.24-30, 37-39, 42
Filter funnels, vacuum	
Borosilicate glass	.80-88
Plastic	.90-91
Polysulfone	.90-91
Stainless steel	.89
Filter holder selection guides	.78-79, 98
Filter holders	
Glass	.80-88
Plastic	.90-91, 102-104
Pressure type	.98-114, 117-118

Stainless steel89, 99-101, 105-114
 Vacuum type80-91
 Filter manifolds92-95
 Filter
 Depth25-30
 Membrane2-16, 19-20
 Qualitative32
 Quantitative33
 Filterceps21
 Filtering flasks86
 Flow cells117-118
 Flow rate graphs121
 Food analysis37-38

G

Gas line holders100
 Gauge, pressure116
 Glass fiber filters26, 28-29
 Glass filter holder accessories86
 Glass filter holders80-88
 Green filters4, 6, 20
 Gridded filters6, 19-20, 40

H

High Pressure filter holder101
 HPLC filters10, 13, 15
 Hybridization membrane35
 Hydrolab manifolds94-95
 Hydrophilic PTFE membrane10, 15-16
 Hydrophobic edge membrane6, 40
 Hydrophobic PTFE membrane11-12, 15-16

I

Indicator paper44
 Individually wrapped filters19-20
 Inline filter holders100-106

L

Litmus paper44
 Low ash filter paper33

M

Manifolds
 Hydrolab94-95
 PVC92-93
 Stainless steel94-95
 Sterility test94-95
 MCE see Mixed cellulose esters
 Media, microbiology22

Membrane
 Edge hydrophobic6
 Black4, 6, 14, 20
 Cellulose acetate7, 20
 Chemical resistance122-123
 Coated9
 Green4, 6, 20
 Gridded6, 19-20
 Hydrophilic PTFE10
 Hydrophobic PTFE11-12
 Hydrophobic edge6
 Microbiology18-20
 Mixed cellulose esters (MCE)4-6, 19-20
 Nitrocellulose, pure35
 Nomenclature3
 Nylon13
 Polycarbonate14
 PTFE10-12
 Selection guide2-3
 Sterile19-20
 Supported12
 Method 1311, EPA28
 Method 1664, EPA93
 Microanalysis filter holders78-85
 Microbiological media22
 Microfiber glass filters26, 28-29
 Milk sediment disk37
 Mixed Cellulose Esters (MCE)
 Membrane, nonsterile4-6
 Membrane, sterile18-20
 Monitors
 37 mm40-41
 NO₂43
 Multimedia pressure sealing filter holder113
 Multistation manifolds92-95

N

Nitrocellulose, pure35
 See also Mixed cellulose esters
 NO₂ badge43
 Nobuto blood filter43
 Nucleic acid hybridization35
 Nylon
 Disposable syringe filters15-16
 Membrane13

O

O-ring chemical compatibility130
 Opticlear membrane6

P

PCTE membrane	.14
PES	.49, 56-58
Petri dishes	.21
PFA	.104
pH paper	.44
Phase fraction collector	.88
Phase separating filters	.38
Plastic filter holders	.90-91, 102-104
Pleated filters	.55-64, 66-68
Polycarbonate membrane filter	.14
Polyester mesh separators	.30
Polyethersulfone (PES)	.49, 56-58
Polyethylene bag	.22
Polypropylene	
Capsule filter	.51
Cartridge filter	.63-65, 69-70
Housing, disposable unit	.15-16
Polysulfone filter holder	.90-91
Prefilters	.24-30
Pressure dispensing vessels	.111-112, 115-116
Presterilized	
Disposable filter units	.15-16
Filters	.18-20
Media	.22
Petri dishes	.21
Protein binding	.4, 35
PTFE	
Capsule	.50
Cartridge	.59-60
Depth filter	.30
Membrane	.10-12
Thimble	.42
Pump, vacuum pressure	.96
PVC	
Extraction manifold	.93
Manifold, standard cup	.92

Q

Qualitative filter paper	.32, 34
Quantitative filter paper	.33-34
Quartz fiber filters	.27-29

R

Reservoirs	
Filter holder with	.111-112
Pressure dispensing vessel	.115-116
Reverse phase membranes	.38

S

Saccharomyces cerevisiae	.18
Sampling bags	.22
Sanitary stainless filter holder	.105, 109-110
Silicone treated cellulose	.38
Single station manifold	.93-95
Soxhlet extraction	.37, 42
Stainless steel filter holders	.89, 99-101, 105-114
Sterile media	.22
Sterile membrane filters	.18-20
Sterility testing	.84, 94-95
Stirred cells	.117-118
Stoppers	.86
Sugar paper	.37
Support pads	.41
Syringe filter	
Disposable	.15, 16
Reuseable holders	.99, 102-104
System design	
TCCLP, apparatus	.111
TCCLP, glass fiber filter	.28-29

T

Test papers	.44
Thimbles	.42
Toxicity testing	.111
Track etch membrane filter	.14
Triton free membrane	.4
Type 25S filter pad	.41

U

Universal pH paper	.44
--------------------	-----

V

Vacuum manifold	.92-95
Vacuum pressure pump	.96

W

Western blots	.35
Whirl-pak sampling bags	.22
Whole range pH paper	.44
Wide mouth polysulfone filter unit	.90-91
Wide mouth pressure vessel	.115-116

Catalog No.	Page	Catalog No.	Page	Catalog No.	Page
01511012	44	311220	87	351600	83
03CP...	16	311225	86	351620	83
03JP...	16	311230	86	351800	109
03NP...	16	311280	81	352100	85
07020010	44	311290	81	352200	85
07020020	44	311300	81	352600	89
13CP...	16	311320	87	353100	94
13HP...	16	311380	81	353130	94
13JP...	16	311390	81	353300	94
13NP...	16	311400	82	353330	94
25AS...	16	311404	86	357000	105
25CS...	16	311420	87	357200	105
25JP...	16	311425	86	357400	106
25NP...	16	311426	86	357500	106
25NS...	16	311430	86	361120	94
268040...	71	311450	82	361130	94
268050...	72	311480	82	361600	112
26807000	76	311490	82	361700	112
269801...	73	311500	82	363120	94
269802...	74	311520	87	363130	94
2S...	38	311540	82	363320	94
300800	107	311580	82	363330	94
300900	107	311590	82	370S...	41
301000	99	311600	83	373130	93
301200	99	311620	87	373330	93
301500	112	311680	82	37AS...	40
301790	113	311690	82	501000	90-91
301900	107	311700	88	501020	90-91
301990	113	311750	88	501030	90-91
302000	107	311760	88	501050	90-91
302100	109	311800	88	501200	102
302200	109	311850	88	501300	102
302300	112	311860	88	50CP...	16
302400	112	313000	94	50JP...	16
302500	107	313200	94	530100	104
302590	113	313400	92	530200	104
302600	107	313600	92	540100	102
302700	109	314100	84	720140	116
302800	109	314110	82	720160	116
304500	100	315110	82	720340	116
304700	100	316110	82	720360	116
304900	101	341000	107	720540	116
311000	80	341100	107	720560	116
311004	86	341200	107	730840	116
311008	86	341300	107	730860	116
311050	80	341400	107	731040	116
311100	80	341500	107	731060	116
311110	80	351120	94	740140	116
311200	81	351130	94	740160	116

INDEX BY CATALOG NUMBER

Catalog No.	Page	Catalog No.	Page	Catalog No.	Page
740340	116	A045T...	19	J010A...	12
740360	116	A045U...	6	J020A...	12
740540	116	A045V...	6	J050A...	12
740560	116	A045W...	20	J100A...	12
740840	116	A065A...	5	K010A...	14
740860	116	A065F...	19	K020A...	14
741040	116	A065G...	19	K020N...	14
741060	116	A065H...	19	K040A...	14
800000	21	A080A...	5	K040N...	14
800100	21	A080B...	6	K080A...	14
800101	21	A080F...	19	K100A...	14
800500	21	A080G...	19	K300A...	14
800501	21	A080H...	19	K800A...	14
800600	43	A080M...	20	MFC200	22
800700	43	A080N...	6	MFCR...	22
800701	43	A080P...	6	MHPC200	22
800800	96	A080Q...	20	MPA200	22
800801	96	A080R...	20	N010A...	13
84...	42	A080V...	6	N022A...	13
86R...	42	A080X...	6, 41	N045A...	13
88R...	42	A100A...	5	N065A...	13
88RH...	42	A100H...	19	N080A...	13
89...	42	A300A...	5	N120A...	13
89S...	42	A500A...	5	N500A...	13
A010A...	5	B200G047A	21	No. 1...	32
A020A...	5	C020A...	8	No. 101...	32
A020C...	19	C020G...	20	No. 1026	37
A020F...	19	C045A...	8	No. 1026-B	37
A020G...	19	C045G...	20	No. 106...	37
A020J...	6	C080A...	8	No. 131...	32
A020K...	6	C300A...	8	No. 1514A...	36
A020S...	19	CCF...	50	No. 151B...	36
A030A...	5	CCFH...	50	No. 1640	25
A045A...	5	CCG...	52	No. 1650...	25
A045B...	6	CCP...	51	No. 2...	32
A045C...	19	CCS...	49	No. 231...	32
A045D...	19	COLI...	22	No. 232...	32
A045E...	19	DP...	26	No. 235...	32
A045F...	19	DS...	30	No. 27...	25
A045G...	19	GA...	26	No. 3...	33
A045H...	19	GB...	26	No. 4A...	33
A045J...	6	GC...	26	No. 50...	36
A045K...	6	GD...	26	No. 51A...	36
A045M...	20	GF...	26	No. 526...	36
A045N...	6	GS...	26	No. 5A	33
A045P...	6	H010A...	10	No. 5B...	33
A045Q...	19, 20	H020A...	10	No. 5C...	33
A045R...	19, 20	H050A...	10	No. 6...	33
A045S...	19	H100A...	10	No. 7...	33

Catalog No.	Page
OH...	69
PF...	30, 38
PG60...	39
QR...	27
S020A...	35
S045A...	35
T010A...	11
T020A...	11
T050A...	11
T080A...	11
T100A...	11
T300A...	11
TC...	66-67
TCF...	59
TCFH...	60
TCG...	68
TCI200	22
TCM200	22
TCP...	63-64
TCPD...	65
TCPE...	64
TCR...	55
TCS...	56
TCS-E...	58
TCS-G...	57
TCY...	61
TCYE...	62
TYPE25S...	41
UHP-13C	46
UNIV...	44
USY...	45
W...	70
WP006ZWS	22
WP018ZWS	22
WP027ZWS	22
WP100MWT	22
WP1254WSP	22
W-R	44
Y008A...	9
Y020A...	9
Y100A...	9
YMM200	22